BOND-PRIME Hematoxylin

Catalog No: AR0096

Leica Biosystems Newcastle Ltd
Balliol Business Park West
Benton Lane
Newcastle Upon Tyne NE12 8EW
United Kingdom
☎ +44 191 215 4242

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Intended Use

For *in vitro* diagnostic use.

BOND-PRIME Hematoxylin is a ready to use reagent that facilitates cell nuclei (blue) counterstaining of formalin-fixed, paraffin-embedded tissue in immunohistochemistry (IHC) and *in situ* hybridization (ISH) using the automated BOND-PRIME staining system.

The clinical interpretation of any staining or its absence should be complemented by morphological studies and proper controls and should be evaluated within the context of the patient's clinical history and other diagnostic tests by a qualified pathologist.

BOND-PRIME Hematoxylin is intended to be used with other devices for visualization of staining and as such, qualitative or semi-quantitative diagnostic function including specific disease indication and intended use population is described within the associated device labelling as required for that device.

Test Principle

Immunohistochemical techniques can be used to demonstrate the presence of antigens in tissue and cells.

In IHC procedures, the Hematoxylin counterstain generates a blue-based stain providing a suitable contrast against non-blue chromogens, aiding in the visualization and localization of the biomarker. Hematoxylin stains cell nuclei and its most common application is the Hematoxylin and Eosin stain.¹

Using BOND-PRIME Hematoxylin in combination with the automated BOND-PRIME system reduces the possibility of human error and inherent variability resulting from individual reagent dilution, manual pipetting and reagent application.

Reagent Description

BOND-PRIME Hematoxylin counterstain contains a <0.1% Hematoxylin solution.

Total volume = 30 mL.

Dilution and Mixing

BOND-PRIME Hematoxylin is ready-to-use. Do not dilute this reagent.

Materials Required But Not Provided

Refer to BOND-PRIME User Manual for a complete list of materials required for specimen treatment and immunohistochemical staining using the BOND-PRIME system.

Storage and Stability

Store at 2–8 °C. Do not use after the expiration date indicated on the container label. Return to 2–8 °C when not in use.

Ensure there are no obvious signs indicating contamination and/or instability of BOND-PRIME Hematoxylin counterstain.

If unexpected staining is observed that cannot be explained by variations in laboratory procedures, and a problem with the product is suspected, contact your local distributor or the regional office of Leica Biosystems immediately.

Storage conditions other than those specified above must be verified by the user.²

Instructions for Use

The BOND-PRIME Hematoxylin was developed for use with the BOND-PRIME Polymer DAB Detection System (DS9284) and BOND-PRIME Polymer AP Detection System (DS9304) on the automated BOND-PRIME system and is the default counterstain reagent for both detection systems.

Precautions and Product Specific Limitations

- To obtain a copy of the Material Safety Data Sheet contact your local distributor or regional office of Leica Biosystems, or alternatively, visit the Leica Biosystems' web site, LeicaBiosystems.com
- Specimens, before and after fixation, and all materials exposed to them, should be handled as if capable of transmitting infection and disposed of with proper precautions.³ Never pipette reagents by mouth and avoid contacting the skin and mucous membranes with reagents or specimens. If reagents or specimens come in contact with sensitive areas, wash with copious amounts of water. Seek medical advice.
- Consult Federal, State or local regulations for disposal of any potentially toxic components.
- Minimize microbial contamination of reagents or an increase in non-specific staining may occur.
- Should any serious incident occur in relation to the product, the user shall report the incident to the manufacturer and the competent authority of the Member State in which the user is established.

Troubleshooting

Refer to reference 1 for remedial action.

Contact your local distributor or the regional office of Leica Biosystems to report unusual staining.

Further Information

Further information on immunostaining with BOND-PRIME reagents can be found in the BOND-PRIME User Manual.

Bibliography

1. Bancroft JD and Stevens A. Theory and Practice of Histological Techniques. 4th Edition. Churchill Livingstone, New York. 1996.
2. Clinical Laboratory Improvement Amendments of 1988, Final Rule 57 FR 7163 February 28, 1992.
3. Villanova PA. National Committee for Clinical Laboratory Standards (NCCLS). Protection of laboratory workers from infectious diseases transmitted by blood and tissue; proposed guideline. 1991; 7(9). Order code M29-P.

Intellectual Property

Copyright © 2020 Leica Biosystems Newcastle Ltd. All rights reserved. LEICA and the Leica Logo are registered trademarks of Leica Microsystems IR GmbH in the US and many other countries. BOND and Novocastra are registered trademarks of the Leica Biosystems group of companies in the USA and optionally in other countries.

Change History

Revision: Date of Issue	Detail of Revision
17 February 2022	<p>Intended Purpose: Update of device Intended Purpose, in accordance with REGULATION (EU) 2017/746 Chapter III 20.4.1.</p> <p>Test Principle: Update of section header from "Principle of Procedure", and removal of 'Leica Biosystems Advanced Staining' term (rewording of section).</p> <p>Reagent Description: Update of section header from "Reagent Composition".</p> <p>Precautions and Product Specific Limitations: Inclusion of statement indicating serious incident protocol.</p> <p>Instructions for Use: Movement of section to above "Precautions and Product Specific Limitations" section.</p> <p>Intellectual Property: Addition of new "Intellectual Property" section.</p> <p>Change History: Addition of new "Change History" section.</p>

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Bibliographie

- 1. Bancroft JD and Stevens A. Theory and Practice of Histological Techniques. 4th Edition. Churchill Livingstone, New York. 1996.
- 2. Clinical Laboratory Improvement Amendments of 1988, Final Rule 57 FR 7163 February 28, 1992.
- 3. Villanova PA. National Committee for Clinical Laboratory Standards (NCCLS). Protection of laboratory workers from infectious diseases transmitted by blood and tissue; proposed guideline. 1991; 7(9). Order code M29-P.

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Copyright© 2020 Leica Biosystems Newcastle, Inc. Tous droits réservés. LEICA et le logo Leica sont des marques déposées de Leica Microsystems IR GmbH aux États-Unis et sont utilisées dans plusieurs autres pays. BOND et Novocastra sont des marques déposées du groupe de sociétés Leica Biosystems aux États-Unis et éventuellement dans d'autres pays.
Copyright © 2020 Leica Biosystems Newcastle Ltd. Tutti i diritti riservati. LEICA e il logo Leica sono marchi registrati di Leica Microsystems IR GmbH negli USA e in diversi altri Paesi. BOND e Novocastra sono marchi registrati del gruppo di aziende Leica Biosystems negli USA e facoltativamente in altri Paesi.
Copyright © 2020 Leica Biosystems Newcastle Ltd. Alle Rechte vorbehalten. LEICA und das Leica-Logo sind eingetragene Marken von Leica Microsystems IR GmbH in den USA und vielen anderen Ländern. BOND und Novocastra sind eingetragene Marken der Leica Biosystems Unternehmensgruppen in den USA und wahlweise in anderen Ländern.
Copyright © 2020 Leica Biosystems Newcastle Ltd. Todos los derechos reservados. LEICA y el logotipo de Leica son marcas comerciales registradas de Leica Microsystems IR GmbH en los Estados Unidos y en muchos otros países. BOND y Novocastra son marcas registradas del grupo de empresas Leica Biosystems en los Estados Unidos y, opcionalmente, en otros países.
Copyright © 2020 Leica Biosystems Imaging, Inc. Todos os direitos reservados. Leica e o logótipo Leica são marcas comerciais registadas da Leica Microsystems IR GmbH nos EUA e noutros países. BOND e Novocastra são marcas comerciais registadas do grupo de empresas Leica Biosystems nos EUA e opcionalmente noutros países.
Copyright © 2020 Leica Biosystems Newcastle Ltd. Med ensamrätt. LEICA och Leica-logotypen är registrerade varumärken som tillhör Leica Microsystems IR GmbH i USA och andra länder. BOND och Novocastra är registrerade varumärken som tillhör företagskoncernen Leica Biosystems i USA och eventuellt i andra länder.