BOND-PRIME Epitope Retrieval Solution 1

Catalog No: AR0086

Leica Biosystems Newcastle Ltd
Balliol Business Park West
Benton Lane
Newcastle Upon Tyne NE12 8EW
United Kingdom
☎ +44 191 215 4242

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Intended Purpose

For *in vitro* diagnostic use.

BOND-PRIME Epitope Retrieval Solution 1 is a ready to use epitope retrieval solution intended for the heat-induced epitope retrieval (HIER) of formalin-fixed, paraffin-embedded tissue using the automated BOND-PRIME system.

The clinical interpretation of any staining or its absence should be complemented by morphological studies using proper controls and should be evaluated within the context of the patient's clinical history and other diagnostic tests by a qualified pathologist.

BOND-PRIME Epitope Retrieval Solution 1 is intended to be used with other devices during immunohistochemistry (IHC) and *in situ* hybridization (ISH) procedures and as such, qualitative or semi-quantitative diagnostic function including specific disease indication and intended use population is described within the associated device labelling as required for that device.

Test Principle

Immunohistochemical techniques can be used to demonstrate the presence of antigens in tissue and cells.

The use of a HIER pretreatment on formalin-fixed, paraffin-embedded tissue on the automated BOND-PRIME system restores epitopes that have been modified by formalin fixation, allowing accessibility of the primary antibody to the epitope.

Reagent Description

BOND-PRIME Epitope Retrieval Solution 1 contains a citrate based buffer and surfactant. Total volume = 1 L.

pH 5.9–6.1 @ 25 °C, as measured at the point of formulation and product release.

Dilution and Mixing

BOND-PRIME Epitope Retrieval Solution 1 is ready to use. Do not dilute this reagent.

Materials Required But Not Provided

Refer to BOND-PRIME User Manual for a complete list of materials required for specimen treatment and immunohistochemical staining using the automated BOND-PRIME system.

Storage and Stability

Store at 2–8 °C prior to pouring into the bulk container on board the BOND-PRIME processing module. Do not use after the expiration date indicated on the bottle label.

The signs indicating contamination and/or instability of BOND-PRIME Epitope Retrieval Solution 1 are turbidity of the solution, odor development, and presence of precipitate.

Storage conditions other than those specified above must be verified by the user.¹

Instructions for Use

Ensure tamper proof seal is intact prior to opening.

The barcode on the bulk bottle of BOND-PRIME Epitope Retrieval Solution 1 should be read by the barcode scanner on the front of the BOND-PRIME Processing Module.

The BOND-PRIME Epitope Retrieval Solution 1 can then be poured directly into the bulk container marked "Epitope Retrieval 1" located within the BOND-PRIME Processing Module.

This container can hold up to 1.25 L.

PLEASE ENSURE THE TAMPER EVIDENT SEAL IS DISPOSED OF IMMEDIATELY AFTER OPENING.

Precautions and Product Specific Limitations

- To obtain a copy of the Material Safety Data Sheet contact your local distributor or regional office of Leica Biosystems, or alternatively, visit the Leica Biosystems' web site, LeicaBiosystems.com
- Specimens, before and after fixation, and all materials exposed to them, should be handled as if capable of transmitting infection and disposed of with proper precautions.² Never pipette reagents by mouth and avoid contacting the skin and mucous membranes with reagents or specimens. If reagents or specimens come in contact with sensitive areas, wash with copious amounts of water. Seek medical advice.
- Consult Federal, State or local regulations for disposal of any potentially toxic components.
- Minimize microbial contamination of reagents or an increase in non-specific staining may occur.
- Should any serious incident occur in relation to the product, the user shall report the incident to the manufacturer and the competent authority of the Member State in which the user is established.

Troubleshooting

Reference 3 may aid in remedial action.

Contact your local distributor or the regional office of Leica Biosystems to report unusual staining.

Further Information

Further information on immunostaining with BOND-PRIME reagents can be found in the BOND-PRIME User Manual.

Bibliography

1. Clinical Laboratory Improvement Amendments of 1988, Final Rule 57 FR 7163 February 28, 1992.
2. Villanova PA. National Committee for Clinical Laboratory Standards (NCCLS). Protection of laboratory workers from infectious diseases transmitted by blood and tissue; proposed guideline. 1991; 7(9). Order code M29-P.
3. Bancroft JD and Stevens A. Theory and Practice of Histological Techniques. 4th Edition. Churchill Livingstone, New York. 1996.

Intellectual Property

Copyright © 2020 Leica Biosystems Newcastle Ltd. All rights reserved. LEICA and the Leica Logo are registered trademarks of Leica Microsystems IR GmbH in the US and many other countries. BOND and Novocastra are registered trademarks of the Leica Biosystems group of companies in the USA and optionally in other countries.

Change History

| Revision: Date of Issue | Detail of Revision |
|-------------------------|---|
| 17 February 2022 | <p>Intended Purpose: Update of device Intended Purpose, in accordance with REGULATION (EU) 2017/746 Chapter III 20.4.1.</p> <p>Test Principle: Update of section header from "Principle of Procedure", and removal of 'Leica Biosystems Advanced Staining' term (rewording of section).</p> <p>Reagent Description: Update of section header from "Reagent Composition".</p> <p>Precautions and Product Specific Limitations: Inclusion of statement indicating serious incident protocol.</p> <p>Instructions for Use: Movement of section to above "Precautions and Product Specific Limitations" section.</p> <p>Intellectual Property: Addition of new "Intellectual Property" section.</p> <p>Change History: Addition of new "Change History" section.</p> |

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- 1. Clinical Laboratory Improvement Amendments of 1988, Final Rule 57 FR 7163 February 28, 1992.
- 2. Villanova PA. National Committee for Clinical Laboratory Standards (NCCLS). Protection of laboratory workers from infectious diseases transmitted by blood and tissue; proposed guideline. 1991; 7(9). Order code M29-P.
- 3. Bancroft JD and Stevens A. Theory and Practice of Histological Techniques. 4th Edition. Churchill Livingstone, New York. 1996.

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- Clinical Laboratory Improvement Amendments of 1988, Final Rule 57 FR 7163 February 28, 1992.
- Villanova PA. National Committee for Clinical Laboratory Standards (NCCLS). Protection of laboratory workers from infectious diseases transmitted by blood and tissue; proposed guideline. 1991; 7(9). Order code M29-P.
- Bancroft JD and Stevens A. Theory and Practice of Histological Techniques. 4th Edition. Churchill Livingstone, New York. 1996.

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Copyright© 2020 Leica Biosystems Newcastle, Inc. Tous droits réservés. LEICA et le logo Leica sont des marques déposées de Leica Microsystems IR GmbH aux États-Unis et sont utilisées dans plusieurs autres pays. BOND et Novocastra sont des marques déposées du groupe de sociétés Leica Biosystems aux États-Unis et éventuellement dans d'autres pays.
Copyright © 2020 Leica Biosystems Newcastle Ltd. Tutti i diritti riservati. LEICA e il logo Leica sono marchi registrati di Leica Microsystems IR GmbH negli USA e in diversi altri Paesi. BOND e Novocastra sono marchi registrati del gruppo di aziende Leica Biosystems negli USA e facoltativamente in altri Paesi.
Copyright © 2020 Leica Biosystems Newcastle Ltd. Alle Rechte vorbehalten. LEICA und das Leica-Logo sind eingetragene Marken von Leica Microsystems IR GmbH in den USA und vielen anderen Ländern. BOND und Novocastra sind eingetragene Marken der Leica Biosystems Unternehmensgruppen in den USA und wahlweise in anderen Ländern.
Copyright © 2020 Leica Biosystems Newcastle Ltd. Todos los derechos reservados. LEICA y el logotipo de Leica son marcas comerciales registradas de Leica Microsystems IR GmbH en los Estados Unidos y en muchos otros países. BOND y Novocastra son marcas registradas del grupo de empresas Leica Biosystems en los Estados Unidos y, opcionalmente, en otros países.
Copyright © 2020 Leica Biosystems Imaging, Inc. Todos os direitos reservados. Leica e o logótipo Leica são marcas comerciais registadas da Leica Microsystems IR GmbH nos EUA e noutros países. BOND e Novocastra são marcas comerciais registadas do grupo de empresas Leica Biosystems nos EUA e opcionalmente noutros países.
Copyright © 2020 Leica Biosystems Newcastle Ltd. Med ensamrätt. LEICA och Leica-logotypen är registrerade varumärken som tillhör Leica Microsystems IR GmbH i USA och andra länder. BOND och Novocastra är registrerade varumärken som tillhör företagskoncernen Leica Biosystems i USA och eventuellt i andra länder.