# BOND-PRIME

## Polymer AP Detection System

**Catalog No: DS9304**

Leica Biosystems Newcastle Ltd  
Balliol Business Park West  
Benton Lane  
Newcastle Upon Tyne NE12 8EW  
United Kingdom  
☎ +44 191 215 4242

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# BOND-PRIME Polymer AP Detection System

## Catalog No: DS9304

For Professional Use Only

### Intended Purpose

For *in vitro* diagnostic use.

BOND-PRIME Polymer AP Detection System is a biotin-free, polymeric alkaline phosphatase (AP)-linker antibody conjugate system intended for the visualization of tissue-bound mouse IgG, mouse IgM and rabbit primary antibodies. It is intended for target visualization by immunohistochemistry (IHC) in sections of formalin fixed, paraffin-embedded tissue using the automated BOND-PRIME system.

The clinical interpretation of any staining or its absence should be complemented by morphological studies and proper controls should be evaluated within the context of the patient's clinical history and other diagnostic tests by a qualified pathologist.

BOND-PRIME Polymer AP Detection System is intended to be used with other devices for visualization of staining and as such, qualitative or semi-quantitative diagnostic function including specific disease indication and intended use population is described within the associated device labelling as required for that device.

### Summary and Explanation

Immunohistochemical (IHC) staining techniques allow for the visualization of antigens via the sequential application of a specific antibody to the antigen (primary antibody), a secondary antibody to the primary antibody and an enzyme complex with a chromogenic substrate with interposed washing steps.

BOND-PRIME Polymer AP Detection System utilizes a novel controlled polymerization technology to prepare polymeric AP-linker antibody conjugates. The detection system avoids the use of streptavidin and biotin, and therefore eliminates non-specific staining as a result of endogenous biotin.

BOND-PRIME Polymer AP Detection System works as follows:

- A user-supplied specific primary antibody is applied.
- Post Primary IgG linker reagent localizes mouse antibodies.
- Poly-AP IgG reagent localizes rabbit antibodies.
- The substrate chromogen, Fast Red, visualizes the complex via a red precipitate.

Using the BOND-PRIME Polymer AP Detection System in combination with the automated BOND-PRIME system reduces the potential for human error and inherent variability resulting from individual reagent dilution, manual pipetting and reagent application.

### Reagents Provided

The reagents provided are sufficient for 100 tests.

1. Post Primary AP (15 mL) Rabbit anti-mouse IgG (<10 µg/mL) in Tris-buffered saline containing 10% (v/v) animal serum and 0.1% ProClin™ 950.
2. Polymer AP (15 mL) Poly-AP anti-rabbit IgG (<25 µg/mL) in Tris-buffered saline containing 10% (v/v) animal serum and 0.09% ProClin™ 950.
3. Red Part A (4.5 mL) Activator containing 0.5% ProClin™ 950.
4. Red Part A (4.5 mL) Activator containing 0.5% ProClin™ 950.
5. Red Part B (1.0 mL) Substrate.
6. Red Part C (1.0 mL) Substrate.
7. Red Part D (32 mL) Buffer solution containing 0.5% ProClin™ 950.
8. Red Part D (32 mL) Buffer solution containing 0.5% ProClin™ 950.

### Dilution and Mixing

BOND-PRIME Polymer AP Detection System is optimized for use on the BOND-PRIME system.

Reconstitution, mixing, dilution, or titration of these reagents is not required.

### Materials Required But Not Provided

BOND-PRIME Dewax Solution (Catalog No. AR0084)

BOND-PRIME Epitope Retrieval Solution (Catalog No. AR0086 or AR0087)

BOND Enzyme Pretreatment Kit (AR9551)

BOND-PRIME Wash Solution Concentrate (Catalog No. AR0085)

BOND-PRIME Hematoxylin (Catalog No. AR0096)

BOND-PRIME system

Refer to BOND-PRIME User Manual for a complete list of materials required for specimen treatment and immunohistochemical staining using the BOND-PRIME system.

### Storage and Stability

Store at 2–8 °C. Do not freeze. Do not use after the expiration date indicated on the tray handle label. Return to 2–8 °C immediately after use.

There are no obvious signs to indicate instability of this product, therefore positive and negative controls should be run simultaneously with unknown specimens (refer to BOND-PRIME User Manual).

If unexpected staining is observed that cannot be explained by variations in laboratory procedures, and a problem with the detection system is suspected, contact your local distributor or the regional office of Leica Biosystems immediately.

Storage conditions other than those specified above must be verified by the user.<sup>1</sup>

## Instructions for Use

The BOND-PRIME Polymer AP Detection System was developed for use on the automated BOND-PRIME system using \*IHC Protocol J. Please refer to individual primary antibody IFUs for default protocol to use with the detection system. Operating parameters for application of the detection system reagents on the BOND-PRIME Processing Module have been optimized by Leica Biosystems. These can be displayed by following the instructions in your BOND-PRIME User Manual.

## Precautions

- Restricted to professional users.
- This detection system is for *in vitro* diagnostic use.
- WARNING - some mounting media may cause the red chromogen to crystallize or fade.
- DO NOT USE DPX, Entellan®, Leica CV Mount, Eukitt® quick-hardening mounting medium, Surgipath MM 24®, Surgipath Sub-X® and Pertex mounting media as they are not suitable for use with this product.
- Fading or crystallization has not been observed with: Leica CV Ultra Mounting Media (Cat 14070937891), Poly-Mount® (Polysciences Inc. Cat 08381), VectaMount™ (Vector Laboratories Cat H-5000) and Limonene Mount (Electron Microscopy Sciences Cat 17987-01).
- The activator and substrate solutions may cause irritation to the skin, eyes, mucous membranes and upper respiratory tract. Wear disposable gloves when handling reagents. Dispose of according to local regulations.
- To obtain a copy of the Material Safety Data Sheet contact your local distributor or regional office of Leica Biosystems. Alternatively, visit the Leica Biosystems' web site, [LeicaBiosystems.com](http://LeicaBiosystems.com)
- Specimens, before and after fixation, and all materials exposed to them, should be handled as if capable of transmitting infection and disposed of with proper precautions.<sup>2</sup> Never pipette reagents by mouth and avoid contacting the skin and mucous membranes with reagents or specimens. If reagents or specimens come in contact with sensitive areas, wash with copious amounts of water. Seek medical advice.
- Consult Federal, State or local regulations for disposal of any potentially toxic components.
- Minimize microbial contamination of reagents or an increase in non-specific staining may occur.
- Incubation times or temperatures other than those specified may give erroneous results. Any such change must be validated by the user.<sup>1</sup>
- Do not mix reagents from different detection systems.
- Should any serious incident occur in relation to the product, the user shall report the incident to the manufacturer and the competent authority of the Member State in which the user is established.

## Product Specific Limitations

BOND-PRIME Polymer AP Detection System has been optimized at Leica Biosystems for use with BOND-PRIME ancillary reagents. Laboratories may use their own primary antibodies provided they have been diluted to an appropriate concentration with BOND Primary Antibody Diluent (Catalog No. AR9352). Users who deviate from recommended test procedures must accept responsibility for interpretation of patient results under these circumstances.

The appropriate concentration of user's own primary antibodies may vary, due to variation in tissue fixation and the effectiveness of antigen enhancement, and must be determined empirically. Negative reagent controls should be used when optimizing retrieval conditions and primary antibody concentrations.

The clinical interpretation of any staining or its absence should be complemented by morphological studies and proper controls. They should be evaluated within the context of the patient's clinical history and other diagnostic tests by a qualified pathologist.

The BOND-PRIME Polymer AP Detection System must be used with laboratory best practice in the use of tissue controls. For assurance, laboratories should stain each patient sample in conjunction with positive, negative and other tissue specific controls as needed.

Some mounting media are not compatible with BOND-PRIME Polymer AP Detection System (see Precautions).

Rapid alcohol/xylene dehydration is recommended.

Contact your local distributor or regional office of Leica Biosystems for further information.

## Troubleshooting

Refer to reference 3 for remedial action.

Contact your local distributor or the regional office of Leica Biosystems to report unusual staining.

## Further Information

Further information on immunostaining with BOND-PRIME reagents can be found in the BOND-PRIME User Manual.

## Bibliography

1. Clinical Laboratory Improvement Amendments of 1988, Final Rule 57 FR 7163 February 28, 1992.
2. Villanova PA. National Committee for Clinical Laboratory Standards (NCCLS). Protection of laboratory workers from infectious diseases transmitted by blood and tissue; proposed guideline. 1991; 7(9). Order Code M29-P.
3. Bancroft JD and Stevens A. Theory and Practice of Histological Techniques. 4th Edition. Churchill Livingstone, New York. 1996.

## Intellectual Property

Copyright © 2020 Leica Biosystems Newcastle Ltd. All rights reserved. LEICA and the Leica Logo are registered trademarks of Leica Microsystems IR GmbH in the US and many other countries. BOND and Novocastra are registered trademarks of the Leica Biosystems group of companies in the USA and optionally in other countries.

Change History

Revision: Date of Issue	Detail of Revision
17 February 2022	<p><b>Intended Purpose:</b> Update of device Intended Purpose, in accordance with REGULATION (EU) 2017/746 Chapter III 20.4.1.</p> <p><b>Summary and Explanation:</b> Update of initial primary usage statement.</p> <p><b>Materials Required But Not Provided:</b> Update to list materials.</p> <p><b>Precautions:</b> Inclusion of statement indicating serious incident protocol.</p> <p><b>Instructions for Use:</b> Movement of section to above "Precautions" section.</p> <p><b>Intellectual Property:</b> Addition of new "Intellectual Property" section.</p> <p><b>Change History:</b> Addition of new "Change History" section.</p>

Revision / Date of Issue

17 February 2022

# BOND-PRIME Polymer AP Detection System

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**Bibliographie**

- 1. Clinical Laboratory Improvement Amendments of 1988, Final Rule 57 FR 7163 February 28, 1992.
- 2. Villanova PA. National Committee for Clinical Laboratory Standards (NCCLS). Protection of laboratory workers from infectious diseases transmitted by blood and tissue; proposed guideline. 1991; 7(9). Order Code M29-P.
- 3. Bancroft JD and Stevens A. Theory and Practice of Histological Techniques. 4th Edition. Churchill Livingstone, New York. 1996.

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## Bibliografi

1. Clinical Laboratory Improvement Amendments of 1988, Final Rule 57 FR 1763 February 28, 1992.
2. Villanova PA. National Committee for Clinical Laboratory Standards (NCCLS). Protection of laboratory workers from infectious diseases transmitted by blood and tissue; proposed guideline. 1991; 7(9). Order Code M29-P.
3. Bancroft JD and Stevens A. Theory and Practice of Histological Techniques. 4th Edition. Churchill Livingstone, New York. 1996.

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Copyright© 2020 Leica Biosystems Newcastle, Inc. Tous droits réservés. LEICA et le logo Leica sont des marques déposées de Leica Microsystems IR GmbH aux États-Unis et sont utilisées dans plusieurs autres pays. BOND et Novocastra sont des marques déposées du groupe de sociétés Leica Biosystems aux États-Unis et éventuellement dans d'autres pays.
Copyright © 2020 Leica Biosystems Newcastle Ltd. Tutti i diritti riservati. LEICA e il logo Leica sono marchi registrati di Leica Microsystems IR GmbH negli USA e in diversi altri Paesi. BOND e Novocastra sono marchi registrati del gruppo di aziende Leica Biosystems negli USA e facoltativamente in altri Paesi.
Copyright © 2020 Leica Biosystems Newcastle Ltd. Alle Rechte vorbehalten. LEICA und das Leica-Logo sind eingetragene Marken von Leica Microsystems IR GmbH in den USA und vielen anderen Ländern. BOND und Novocastra sind eingetragene Marken der Leica Biosystems Unternehmensgruppen in den USA und wahlweise in anderen Ländern.
Copyright © 2020 Leica Biosystems Newcastle Ltd. Todos los derechos reservados. LEICA y el logotipo de Leica son marcas comerciales registradas de Leica Microsystems IR GmbH en los Estados Unidos y en muchos otros países. BOND y Novocastra son marcas registradas del grupo de empresas Leica Biosystems en los Estados Unidos y, opcionalmente, en otros países.
Copyright © 2020 Leica Biosystems Imaging, Inc. Todos os direitos reservados. Leica e o logótipo Leica são marcas comerciais registadas da Leica Microsystems IR GmbH nos EUA e noutros países. BOND e Novocastra são marcas comerciais registadas do grupo de empresas Leica Biosystems nos EUA e opcionalmente noutros países.
Copyright © 2020 Leica Biosystems Newcastle Ltd. Med ensamrätt. LEICA och Leica-logotypen är registrerade varumärken som tillhör Leica Microsystems IR GmbH i USA och andra länder. BOND och Novocastra är registrerade varumärken som tillhör företagskoncernen Leica Biosystems i USA och eventuellt i andra länder.